



FACILITIES
Emtech Laboratories
7745 Garland Circle
Roanoke, VA 24019

MAILING
Emtech Laboratories
P.O. Box 12900
Roanoke, VA 24022

PHONE
Local: (540) 265-9156
Toll Free: (800) 336-5719
Fax: (540) 265-9164

Half Shell

This style is a modified shell. The anti-helix and parts of the upper concha of the mold are removed. It provides better retention than the canal mold and is used in many cases when the canal mold will not provide proper retention. Also make it a little easier to insert and remove from the ear.

