



**FACILITIES**  
Emtech Laboratories  
7745 Garland Circle  
Roanoke, VA 24019

**MAILING**  
Emtech Laboratories  
P.O. Box 12900  
Roanoke, VA 24022

**PHONE**  
Local: (540) 265-9156  
Toll Free: (800) 336-5719  
Fax: (540) 265-9164

## Innovative Tube Lock

Innovative tube lock has a unique flanged tube lock designed to work with any silicone ear piece to secure tubing without the need for cement. Emtech pre-attaches the tubing to the lock so you have no mess and no waiting for cement to dry.

Innovative tube locks can be applied to any silicone ear piece or purchased separately as replacement items. Also available in Emdry treated.

