



FACILITIES
Emtech Laboratories
7745 Garland Circle
Roanoke, VA 24019

MAILING
Emtech Laboratories
P.O. Box 12900
Roanoke, VA 24022

PHONE
Local: (540) 265-9156
Toll Free: (800) 336-5719
Fax: (540) 265-9164

Receiver

Ordinarily used with hearing instruments having external receiver, which is snapped into the clip ring of the earpiece. It fills the entire concha and is used mostly for body type instruments. You may choose any material or color for this style. The receiver is also used for auditory trainer molds, pilot molds, and anesthesiologist earpieces.

