



FACILITIES

Emtech Laboratories
7745 Garland Circle
Roanoke, VA 24019

MAILING

Emtech Laboratories
P.O. Box 12900
Roanoke, VA 24022

PHONE

Local: (540) 265-9156
Toll Free: (800) 336-5719
Fax: (540) 265-9164

ER15- Custom Musician Earplugs

Custom Musician Earplugs have 3 choices of interchangeable filters. ER9, ER15 or ER25.

The filters come in the standard style, partially countersunk or completely countersunk. When using our Musician's earplugs sound quality is clearer and more natural and fatigue associated with noise exposure is reduced. The world doesn't sound muffled.

Musician earplugs are also great for those who are exposed to 90-120dB sound levels for various time periods and need to hear accurately such as:

Dentists, Pilots and airline personnel who do not require maximum attenuation, recording engineers, sound crews, most industrial workers, machinists, restaurant and dance club employees and rock music listeners.

